

## 1640 无糖液体培养基（含 L-谷氨酰胺，不含丙酮酸钠、HEPES）

货号：BL1126A

规格：500ml

保存温度：2-8℃，避光保存，有效期 12 个月。

Components	Molecular Weight	Concentration (mg/L)	mM
<b>Amino Acids</b>			
Glycine	75	10	0.13333334
L-Arginine	174	200	1.1494253
L-Asparagine	132	50	0.37878788
L-Aspartic acid	133	20	0.15037593
L-Cystine 2HCl	313	65	0.20766774
L-Glutamic Acid	147	20	0.13605443
L-Glutamine	146	300	2.0547945
L-Histidine	155	15	0.09677419
L-Hydroxyproline	131	20	0.15267175
L-Isoleucine	131	50	0.3816794
L-Leucine	131	50	0.3816794
L-Lysine hydrochloride	146	40	0.2739726
L-Methionine	149	15	0.10067114
L-Phenylalanine	165	15	0.09090909
L-Proline	115	20	0.17391305
L-Serine	105	30	0.2857143
L-Threonine	119	20	0.16806723
L-Tryptophan	204	5	0.024509804
L-Tyrosine disodium salt dihydrate	261	29	0.11111111
L-Valine	117	20	0.17094018
<b>Vitamins</b>			
Biotin	244	0.2	0.000819672
Choline chloride	140	3	0.021428572
D-Calcium pantothenate	477	0.25	0.000524109
Folic Acid	441	1	0.002267574
Niacinamide	122	1	0.008196721
Para-Aminobenzoic Acid	137	1	0.00729927
Pyridoxine hydrochloride	206	1	0.004854369
Riboflavin	376	0.2	0.000531915
Thiamine hydrochloride	337	1	0.002967359
Vitamin B12	1355	0.005	3.69E-06
i-Inositol	180	35	0.19444445
<b>Inorganic Salts</b>			

Note: For in vitro research use only, not for diagnostic or therapeutic use, This product is not a medical device.  
 注意：在体外研究使用，不用于诊断或治疗用途，本产品不是医疗装置。



Calcium nitrate (Ca(NO <sub>3</sub> ) <sub>2</sub> 4H <sub>2</sub> O)	236	100	0.42372882
Magnesium Sulfate (MgSO <sub>4</sub> ) (anhyd.)	120	48.84	0.407
Potassium Chloride (KCl)	75	400	5.3333335
Sodium Bicarbonate (NaHCO <sub>3</sub> )	84	2000	23.809525
Sodium Chloride (NaCl)	58	6000	103.44827
Sodium Phosphate dibasic (Na <sub>2</sub> HPO <sub>4</sub> ·7H <sub>2</sub> O)	268	1512	5.641791
<b>Other Components</b>			
Glutathione (reduced)	307	1	0.003257329
Phenol Red	376.4	5	0.013283741

Note: For in vitro research use only, not for diagnostic or therapeutic use, This product is not a medical device.  
 注意: 在体外研究使用, 不用于诊断或治疗用途, 本产品不是医疗装置。

