

DMEM 无糖液体培养基（不含 L-谷氨酰胺，丙酮酸钠、HEPES、酚红）

货号：BL1125A

规格：500ml

保存温度：2-8℃，避光保存，有效期 12 个月。

Components	Molecular Weight	Concentration (mg/L)	mM
Amino Acids			
Glycine	75	30	0.4
L-Arginine hydrochloride	211	84	0.39810428
L-Cystine 2HCl	313	63	0.20127796
L-Histidine hydrochloride-H ₂ O	210	42	0.2
L-Isoleucine	131	105	0.8015267
L-Leucine	131	105	0.8015267
L-Lysine hydrochloride	183	146	0.7978142
L-Methionine	149	30	0.20134228
L-Phenylalanine	165	66	0.4
L-Serine	105	42	0.4
L-Threonine	119	95	0.79831934
L-Tryptophan	204	16	0.078431375
L-Tyrosine disodium salt dihydrate	261	104	0.39846742
L-Valine	117	94	0.8034188
Vitamins			
Choline chloride	140	4	0.028571429
D-Calcium pantothenate	477	4	0.008385744
Folic Acid	441	4	0.009070295
Niacinamide	122	4	0.032786883
Pyridoxine hydrochloride	204	4	0.019607844
Riboflavin	376	0.4	0.00106383
Thiamine hydrochloride	337	4	0.011869436
i-Inositol	180	7.2	0.04
Inorganic Salts			
Calcium Chloride (CaCl ₂) (anhyd.)	111	200	1.8018018
Ferric Nitrate (Fe(NO ₃) ₃ ·9H ₂ O)	404	0.1	2.48E-04
Magnesium Sulfate (MgSO ₄) (anhyd.)	120	97.67	0.8139166
Potassium Chloride (KCl)	75	400	5.3333335
Sodium Bicarbonate (NaHCO ₃)	84	3700	44.04762
Sodium Chloride (NaCl)	58	6400	110.344826
Sodium Phosphate monobasic (NaH ₂ PO ₄ ·H ₂ O)	138	125	0.9057971

Note: For in vitro research use only, not for diagnostic or therapeutic use, This product is not a medical device.
 注意：在体外研究使用，不用于诊断或治疗用途，本产品不是医疗装置。

